

HUGE

Clear Guide Visible Surgery



Fracture Resistance

- Drill position visible in real time
- Early error detection



High Toughness

- Resists fracture under load
- Avoids brittle failure for more controlled surgery



Impact Resistance

- Resists fracture from drops or force



High Rigidity

- Accurate seating; drill sleeves stay in place

6–8%

Elongation

90–95 MPa

Flexural strength

2200–2400 MPa

Flexural modulus

55–60 MPa

Tensile strength

25–30 J/m

Notched impact

HUGE



Available Options

- Color: Clear
- Pack sizes: 500g / 1000g

Certification

MDR certified

FDA listing

Compatible with mainstream 385/405 nm LCD/DLP printers

Validated printers

- ✓ Asiga Max 2
- ✓ Rayshape Edge Mini
- ✓ Rayshape Edge E2

Validated curing units


- ✓ Creality UV 02
- ✓ Rayshape ShapeCure D

Designed for the HUGE Digital Denture Workflow

- Accurate Fit from Digital Design to Surgery
- Stable Guidance Under Load
- Free from offensive odors or tastes



Accurate Design



Precise Printing



Confident Guidance



Shandong Huge Dental Material Corporation

Address: No. 68 Shanhai Road, Donggang District, Rizhao City,
Shandong Province, 276800, P.R. China. Tel: +86 (633) 2277268
marketing@hugedental.com www.hugedental.com

Rev 10/2026



Follow @Huge Dental